

ETV6 antibody

Product Information

Catalog No.:	FNab02882
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcriptional repressor; binds to the DNA sequence 5'-CCGGAAGT-3'. Plays a role in hematopoiesis and malignant transformation.

Immunogen information

Immunogen:	ets variant 6
Synonyms:	Transcription factor ETV6 ETS translocation variant 6 ETS-related protein Tel1 (Tel) ETV6 TEL TEL1
Observed MW:	57 kDa
Uniprot ID :	P41212

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



A431 cells were subjected to SDS PAGE followed by western blot with FNab02882(ETV6 antibody) at dilution of 1:1000