

EXTL2 antibody

Product Information

Catalog No.:	FNab02911
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

Glycosyltransferase required for the biosynthesis of heparan-sulfate and responsible for the alternating addition of beta-1-4-linked glucuronic acid(GlcA) and alpha-1-4-linked N-acetylglucosamine(GlcNAc) units to nascent heparan sulfate chains.


Immunogen information

Immunogen:	exostoses(multiple)-like 2
Synonyms:	Exostosin-like 2 Alpha-1,4-N-acetylhexosaminyltransferase EXTL2 Alpha-GalNAcT EXTL2 EXT-related protein 2 Glucuronyl- galactosyl-proteoglycan 4-alpha-N- acetylglucosaminyltransferase Processed exostosin-like 2 EXTL2 EXTR2
Observed MW:	37 kDa
Uniprot ID :	Q9UBQ6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:3000

Image:



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02911(EXTL2 antibody) at dilution of 1:400