

F5 antibody

Product Information

Catalog No.:	FNab02923
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Central regulator of hemostasis. It serves as a critical cofactor for the prothrombinase activity of factor Xa that results in the activation of prothrombin to thrombin.

Immunogen information

Immunogen:	coagulation factor V(proaccelerin, labile factor)
Synonyms:	Coagulation factor V Activated protein C cofactor Proaccelerin, labile factor Coagulation factor V heavy chain Coagulation factor V light chain F5
Observed MW:	145 kDa
Uniprot ID :	P12259

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab02923(F5 antibody) at dilution of 1:1000