

PTK2 antibody

Product Information

Catalog No.:	FNab02950
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	PTK2 protein tyrosine kinase 2
Synonyms:	Focal adhesion kinase 1 (FADK 1) Focal adhesion kinase-related nonkinase (FRNK) Protein phosphatase 1 regulatory subunit 71 (PPP1R71) Protein-tyrosine kinase 2 p125FAK pp125FAK PTK2 FAK FAK1
Observed MW:	120 kDa
Uniprot ID :	Q05397

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue using FNab02950(FAK antibody) at dilution of 1:200



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab02950(FAK antibody) at dilution of 1:2000