

RTL8A antibody

Product Information

Catalog No.:	FNab02963
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	family with sequence similarity 127, member B
Synonyms:	Retrotransposon Gag-like protein 8A Mammalian retrotransposon derived protein 8A RTL8A CXX1B FAM127B MAR8A
Observed MW:	13 kDa
Uniprot ID :	Q9BWD3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon using FNab02963(FAM127B antibody) at dilution of 1:100



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab02963(FAM127B antibody) at dilution of 1:500