

HAUS6 antibody

Product Information

Catalog No.:	FNab02980
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	family with sequence similarity 29, member A
Synonyms:	HAUS augmin-like complex subunit 6 HAUS6 DGT6 FAM29A KIAA1574
Observed MW:	109 kDa
Uniprot ID :	Q7Z4H7

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab02980(FAM29A antibody) at dilution of 1:500