

## FAM50A antibody

### Product Information

Catalog No.:	FNab02988
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be a DNA-binding protein or transcriptional factor.

### Immunogen information

Immunogen:	family with sequence similarity 50, member A
Synonyms:	Protein FAM50A Protein HXC-26 Protein XAP-5 FAM50A DXS9928E HXC26 XAP5
Observed MW:	37-40 kDa
Uniprot ID :	Q14320

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02988(FAM50A antibody) at dilution of 1:1500