

## FCN3 antibody

### Product Information

Catalog No.:	FNab03061
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May function in innate immunity through activation of the lectin complement pathway. Calcium-dependent and GlcNAc-binding lectin. Has affinity with GalNAc, GlcNAc, D-fucose, as mono/oligosaccharide and lipopolysaccharides from S.typhimurium and S.minnesota.

### Immunogen information

Immunogen:	ficolin(collagen/fibrinogen domain containing) 3
Synonyms:	Ficolin-3 Collagen/fibrinogen domain-containing lectin 3 p35 Collagen/fibrinogen domain-containing protein 3 Hakata antigen FCN3 FCNH HAKA1
Observed MW:	35 kDa, 50 kDa
Uniprot ID :	O75636

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human placenta tissue were subjected to SDS  
PAGE followed by western blot with  
FNab03061(FCN3 antibody) at dilution of 1:500