

FCRLB antibody

Product Information

Catalog No.:	FNab03062
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Fc receptor-like B
Synonyms:	Fc receptor-like B Fc receptor homolog expressed in B-cells protein 2 (FREB-2) Fc receptor-like and mucin-like protein 2 Fc receptor-like protein 2 Fc receptor-related protein Y (FcRY) FCRLB FCRL2 FCRLM2 FCRY FREB2
Observed MW:	57-60 kDa
Uniprot ID :	Q6BAA4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab03062(FCRLB antibody) at dilution of 1:100



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03062(FCRLB Antibody) at dilution of 1:1000