

FGFBP2 antibody

Product Information

Catalog No.:	FNab03097
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	fibroblast growth factor binding protein 2
Synonyms:	Fibroblast growth factor-binding protein 2 (FGF-BP2, FGF-binding protein 2, FGFBP-2) 37 kDa killer-specific secretory protein (Ksp37) HBp17-related protein (HBp17-RP) FGFBP2 KSP37
Observed MW:	37 kDa
Uniprot ID :	Q9BYJ0

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab03097(FGFBP2 antibody) at dilution of 1:400