

## FREM1 antibody

### Product Information

Catalog No.:	FNab03218
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Extracellular matrix protein that plays a role in epidermal differentiation and is required for epidermal adhesion during embryonic development.

### Immunogen information


Immunogen:	FRAS1 related extracellular matrix 1
Synonyms:	FRAS1-related extracellular matrix protein 1 Protein QBRICK FREM1 C9orf143 C9orf145 C9orf154
Observed MW:	40-43 kDa;70-80 kDa
Uniprot ID :	Q5H8C1

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab03218(FREM1 Antibody) at dilution of 1:50



human plasma tissue were subjected to SDS PAGE followed by western blot with FNab03218(FREM1 Antibody) at dilution of 1:300