

FXYD4 antibody

Product Information

Catalog No.:	FNab03260
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	FXYD domain containing ion transport regulator 4
Synonyms:	FXYD domain-containing ion transport regulator 4 FXYD4
Observed MW:	10 kDa
Uniprot ID :	P59646

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse testis tissue slide using FNa03260(FXYD4 Antibody) at dilution of 1:50