

LGALS2 antibody

Product Information

Catalog No.:	FNab03310
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	lectin, galactoside-binding, soluble, 2
Synonyms:	Galectin-2 (Gal-2) Beta-galactoside-binding lectin L-14-II HL14 Lactose-binding lectin 2 S-Lac lectin 2 LGALS2
Observed MW:	15 kDa
Uniprot ID :	P05162

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab03310(LGALS2 antibody) at dilution of 1:100



human colon tissue were subjected to SDS PAGE followed by western blot with FNab03310(LGALS2 antibody) at dilution of 1:300