

## GFRA3 antibody

### Product Information

Catalog No.:	FNab03440
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Receptor for the glial cell line-derived neurotrophic factor, ARTN(artemin). Mediates the artemin-induced autophosphorylation and activation of the RET receptor tyrosine kinase.

### Immunogen information

Immunogen:	GDNF family receptor alpha 3
Synonyms:	GDNF family receptor alpha-3 (GDNF receptor alpha-3, GDNFR-alpha-3, GFR-alpha-3) GFRA3
Observed MW:	45 kDa
Uniprot ID :	O60609

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab03440( GFRA3 Antibody) at dilution of 1:600