

GIMAP2 antibody

Product Information

Catalog No.:	FNab03456
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


The heterodimer formed by GIMAP2 and GIMAP7 has GTPase activity. In contrast, GIMAP2 has no GTPase activity by itself.

Immunogen information

Immunogen:	GTPase, IMAP family member 2
Synonyms:	GTPase IMAP family member 2 Immunity-associated protein 2 (hIMAP2) GIMAP2 IMAP2
Observed MW:	28 kDa
Uniprot ID :	Q9UG22

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human heart tissue were subjected to SDS PAGE followed by western blot with FNab03456(GIMAP2 antibody) at dilution of 1:300