

GORAB antibody

Product Information

Catalog No.:	FNab03563
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	golgin, RAB6-interacting
Synonyms:	RAB6-interacting golgin N-terminal kinase-like-binding protein 1 (NTKL-BP1, NTKL-binding protein 1, hNTKL-BP1) SCY1-like 1-binding protein 1 (SCYL1-BP1, SCYL1-binding protein 1) GORAB NTKLBP1 SCYL1BP1
Observed MW:	45-50 kDaa
Uniprot ID :	Q5T7V8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Hela cells were subjected to SDS PAGE followed by western blot with FNab03563(GORAB antibody) at dilution of 1:600