

NUP210 antibody

Product Information

Catalog No.:	FNab03571
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	nucleoporin 210kDa
Synonyms:	Nuclear pore membrane glycoprotein 210 (Nuclear pore protein gp210) Nuclear envelope pore membrane protein POM 210 (POM210) Nucleoporin Nup210 Pore membrane protein of 210 kDa NUP210 KIAA0906
Observed MW:	205 kDa
Uniprot ID :	Q8TEM1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab03571(gp210 antibody) at dilution of 1:1000