

GPBP1 antibody

Product Information

Catalog No.:	FNab03577
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Functions as a GC-rich promoter-specific transactivating transcription factor.

Immunogen information


Immunogen:	GC-rich promoter binding protein 1
Synonyms:	Vasculin GC-rich promoter-binding protein 1 Vascular wall-linked protein GPBP1 GPBP SSH6
Observed MW:	45 kDa
Uniprot ID :	Q86WP2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human small intestine slide using FNab03577(GPBP1 Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03577(GPBP1 antibody) at dilution of 1:500