

## GPRASP1-Specific antibody

### Product Information

Catalog No.:	FNab03612
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	G protein-coupled receptor associated sorting protein 1
Synonyms:	G-protein coupled receptor-associated sorting protein 1 (GASP-1) GPRASP1 GASP KIAA0443
Observed MW:	157 kDa
Uniprot ID :	Q5JY77

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab03612(GPRASP1-Specific antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03612(GPRASP1-Specific antibody) at dilution of 1:500