

GRAMD4 antibody

Product Information

Catalog No.:	FNab03629
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role as a mediator of E2F1-induced apoptosis in the absence of p53/TP53.

Immunogen information

Immunogen:	GRAM domain containing 4
Synonyms:	GRAM domain-containing protein 4 Death-inducing protein GRAMD4 DIP KIAA0767
Observed MW:	66, 140 kDa
Uniprot ID :	Q6IC98

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab03629(GRAMD4 antibody) at dilution of 1:1200