

GRN antibody

Product Information

Catalog No.:	FNab03632
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None


Immunogen information

Immunogen:	granulin
Synonyms:	Progranulin (PGRN) Acrogranin Epithelin precursor Glycoprotein of 88 Kda (GP88, Glycoprotein 88) Granulin precursor PC cell-derived growth factor (PCDGF) Proepithelin (PEPI) Paragranulin Granulin-1 Alternative names: Granulin G Granulin-2 Alternative names: Granulin F Granulin-3 Alternative names: Epithelin-2 Granulin B Granulin-4 Alternative names: Epithelin-1 Granulin A Granulin-5 Alternative names: Granulin C Granulin-6 Alternative names: Granulin D Granulin-7 Alternative names: Granulin E GRN
Observed MW:	78 kDa
Uniprot ID :	P28799


Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200

Image:



Immunohistochemistry of paraffin-embedded human heart using FNab03632(PCDGF,GRN Antibody) at dilution of 1:50



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab03632(PCDGF,GRN antibody) at dilution of 1:1000