

HBM-Specific antibody

Product Information

Catalog No.:	FNab03773
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	hemoglobin, mu
Synonyms:	Hemoglobin subunit mu Hemoglobin mu chain Mu-globin HBM HBAP2
Observed MW:	16 kDa
Uniprot ID :	Q6B0K9

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:1000
Image:	



human heart tissue were subjected to SDS PAGE followed by western blot with FNab03773(HBM-Specific antibody) at dilution of 1:500