

## HEMGN antibody

### Product Information

Catalog No.:	FNab03824
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Regulates the proliferation and differentiation of hematopoietic cells. Overexpression block the TPA-induced megakaryocytic differentiation in the K562 cell model. May also prevent cell apoptosis through the activation of the nuclear factor-kappa B(NF-kB).

### Immunogen information

Immunogen:	hemogen
Synonyms:	Hemogen Erythroid differentiation-associated gene protein (EDAG-1) Hemopoietic gene protein Negative differentiation regulator protein HEMGN EDAG NDR
Observed MW:	60 kDa
Uniprot ID :	Q9BXL5

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab03824(HEMGN antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03824(HEMGN antibody) at dilution of 1:500