

---

## HBE1-Specific antibody

### Product Information

Catalog No.:	FNab03827
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	hemoglobin, epsilon 1
Synonyms:	Hemoglobin subunit epsilon Epsilon-globin Hemoglobin epsilon chain HBE1 HBE
Observed MW:	15 kDa
Uniprot ID :	P02100

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



K562 cells were subjected to SDS PAGE followed by western blot with FNab03827(Hemoglobin epsilon-Specific Antibody) at dilution of 1:600