

Histone H2A.z antibody

Product Information

Catalog No.:	FNab03886
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	H2A histone family, member Z
Synonyms:	Histone H2A.Z (H2A/z) H2AZ1 H2AFZ H2AZ
Observed MW:	14 kDa
Uniprot ID :	P0C0S5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab03886(Histone H2A.z antibody) at dilution of 1:300