

HMBOX1 antibody

Product Information

Catalog No.:	FNab03917
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcription factor. Isoform 1 acts as a transcriptional repressor. Isoform 4 has very low activity as a transcriptional repressor.

Immunogen information

Immunogen:	homeobox containing 1
Synonyms:	Homeobox-containing protein 1 Homeobox telomere-binding protein 1 Telomere-associated homeobox-containing protein 1 HMBOX1 HOT1 TAH1
Observed MW:	47 kDa
Uniprot ID :	Q6NT76

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab03917(HMBOX1 antibody) at dilution of 1:1000