

HNRNPUL2 antibody

Product Information

Catalog No.:	FNab03961
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	heterogeneous nuclear ribonucleoprotein U-like 2
Synonyms:	Heterogeneous nuclear ribonucleoprotein U-like protein 2 Scaffold-attachment factor A2 (SAF-A2) HNRNPUL2 HNRPUL2
Observed MW:	110 kDa
Uniprot ID :	Q1KMD3

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:1000; IF: 1:50-1:500
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab03961(HNRNPUL2-Specific Antibody) at dilution of 1:300