

HOMEZ antibody

Product Information

Catalog No.:	FNab03966
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function as a transcriptional regulator.

Immunogen information


Immunogen:	homeobox and leucine zipper encoding
Synonyms:	Homeobox and leucine zipper protein Homez Homeodomain leucine zipper-containing factor HOMEZ KIAA1443
Observed MW:	65 kDa
Uniprot ID :	Q8IX15

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis slide using FNab03966(HOMEZ Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03966(HOMEZ Antibody) at dilution of 1:1000