

HOXA11 antibody

Product Information

Catalog No.:	FNab03973
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Immunogen information

Immunogen:	homeobox A11
Synonyms:	Homeobox protein Hox-A11 Homeobox protein Hox-11 HOXA11 HOX11
Observed MW:	37 kDa
Uniprot ID :	P31270

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03973(HOXA11 antibody) at dilution of 1:1000