

HOXC4 antibody

Product Information

Catalog No.:	FNab03982
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Immunogen information

Immunogen:	homeobox C4
Synonyms:	Homeobox protein Hox-C4 Homeobox protein CP19 Homeobox protein Hox-3E HOXC4 HOX3E
Observed MW:	30 kDa
Uniprot ID :	P09017

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03982(HOXC4 antibody) at dilution of 1:500