

HUS1B antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab04083 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information

| | |
|--------------|---|
| Immunogen: | HUS1 checkpoint homolog b(S. pombe) |
| Synonyms: | Checkpoint protein HUS1B (hHUS1B) HUS1B |
| Observed MW: | 31-35 kDa |
| Uniprot ID : | Q8NHY5 |

Application

| | |
|---------------------------|--------------|
| Reactivity: | Human, Rat |
| Tested Application: | ELISA, WB |
| Recommended dilution: WB: | 1:500-1:2000 |
| Image: | |



RAW264.7 cells were subjected to SDS PAGE followed by western blot with FNab04083(HUS1B antibody) at dilution of 1:500