

## ICP27 antibody

### Product Information

Catalog No.:	FNab04112
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	ICP27
Synonyms:	mRNA export factor Immediate-early protein IE63 Infected cell protein 27 (ICP27) VMW63
Observed MW:	55 kDa
Uniprot ID :	P10238

### Application

Reactivity:	HSV
Tested Application:	ELISA