

IDH1 antibody

Product Information

Catalog No.:	FNab04118
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	isocitrate dehydrogenase 1(NADP+), soluble
Synonyms:	Isocitrate dehydrogenase [NADP] cytoplasmic (IDH, IDH1) Cytosolic NADP-isocitrate dehydrogenase IDPc NADP(+)-specific ICDH Oxalosuccinate decarboxylase IDH1 PICD
Observed MW:	46 kDa
Uniprot ID :	O75874

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using FNab04118(IDH1 Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab04118(IDH1 antibody) at dilution of 1:2000