

IFRD2 antibody

Product Information

Catalog No.:	FNab04152
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	interferon-related developmental regulator 2
Synonyms:	Interferon-related developmental regulator 2 Protein SKMC15 IFRD2
Observed MW:	50 kDa
Uniprot ID :	Q12894

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04152(IFRD2 Antibody) at dilution of 1:600