

## IFT140 antibody

### Product Information

Catalog No.:	FNab04154
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Component of the IFT complex A(IFT-A), a complex required for retrograde ciliary transport. Plays a pivotal role in proper development and function of ciliated cells. Involved in ciliogenesis and cilia maintenance.

### Immunogen information

Immunogen:	intraflagellar transport 140 homolog(Chlamydomonas)
Synonyms:	Intraflagellar transport protein 140 homolog WD and tetratricopeptide repeats protein 2 IFT140 KIAA0590 WDTC2
Observed MW:	140 kDa
Uniprot ID :	Q96RY7

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab04154( IFT140 Antibody) at dilution of 1:300