

IGSF11 antibody

Product Information

Catalog No.:	FNab04190
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Functions as a cell adhesion molecule through homophilic interaction. Stimulates cell growth.

Immunogen information


Immunogen:	immunoglobulin superfamily, member 11
Synonyms:	Immunoglobulin superfamily member 11 (IgSF11) Brain and testis-specific immunoglobulin superfamily protein (Bt-IGSF) V-set and immunoglobulin domain-containing protein 3 IGSF11 BTIGSF CXADRL1 VSIG3
Observed MW:	40-50 kDa
Uniprot ID :	Q5DX21

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab04190(IGSF11 antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab04190(IGSF11 antibody) at dilution of 1:500