

immunoglobulin kappa antibody

Product Information

Catalog No.:	FNab04295
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	immunoglobulin kappa locus
Synonyms:	IGK IGK- IGKC immunoglobulin kappa locus IgG Light Chain (Kappa)
Observed MW:	26 kDa
Uniprot ID :	None

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human spleen using FNab04295(IGK-antibody) at dilution of 1:100



human plasma(0.3ug) tissue were subjected to SDS PAGE followed by western blot with FNab04295(IgG light chain(Kappa) Antibody) at dilution of 1:1000