

KPNB1 antibody

Product Information

Catalog No.:	FNab04305
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	karyopherin(importin) beta 1
Synonyms:	Importin subunit beta-1 Importin-90 Karyopherin subunit beta-1 Nuclear factor p97 Pore targeting complex 97 kDa subunit (PTAC97) KPNB1 NTF97
Observed MW:	95 kDa
Uniprot ID :	Q14974

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:1000;IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab04305(KPNB1 antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab04305(KPNB1 Antibody) at dilution of 1:600