

INHBA antibody

Product Information

Catalog No.:	FNab04318
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	0E10
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	inhibin, beta A
Synonyms:	Inhibin beta A chain Activin beta-A chain Erythroid differentiation protein (EDF) INHBA
Observed MW:	47 kDa
Uniprot ID :	P08476

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab04318 (INHBA antibody) at dilution of 1:200



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab04318 (INHBA antibody) at dilution of 1:1000