

INSRR antibody

Product Information

Catalog No.:	FNab04337
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Receptor with tyrosine-protein kinase activity. Functions as a pH sensing receptor which is activated by increased extracellular pH. Activates an intracellular signaling pathway that involves IRS1 and AKT1/PKB.

Immunogen information

Immunogen:	insulin receptor-related receptor
Synonyms:	Insulin receptor-related protein (IRR) IR-related receptor Insulin receptor-related protein alpha chain Insulin receptor-related protein beta chain INSRR IRR
Observed MW:	150-155 kDa
Uniprot ID :	P14616

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab04337(INSRR antibody) at dilution of 1:500