

ITM2A antibody

Product Information

Catalog No.:	FNab04417
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	integral membrane protein 2A
Synonyms:	Integral membrane protein 2A Protein E25 ITM2A
Observed MW:	45-48 kDa, 30 kDa
Uniprot ID :	O43736

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab04417(ITM2A antibody) at dilution of 1:50