

## KCNJ11 antibody

### Product Information

Catalog No.:	FNab04588
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	potassium inwardly-rectifying channel, subfamily J, member 11
Synonyms:	ATP-sensitive inward rectifier potassium channel 11 IKATP Inward rectifier K(+) channel Kir6.2 Potassium channel, inwardly rectifying subfamily J member 11 KCNJ11
Observed MW:	48 kDa
Uniprot ID :	Q14654

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse spleen tissue were subjected to SDS PAGE followed by western blot with FNab04588(Kir6.2 antibody) at dilution of 1:300