

KLHDC4 antibody

Product Information

Catalog No.:	FNab04607
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	kelch domain containing 4
Synonyms:	Kelch domain-containing protein 4 KLHDC4
Observed MW:	58 kDa
Uniprot ID :	Q8TBB5

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab04607(KLHDC4 antibody) at dilution of 1:500