

KMO antibody

Product Information

Catalog No.:	FNab04668
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	kynurenine 3-monooxygenase(kynurenine 3-hydroxylase)
Synonyms:	Kynurenine 3-monooxygenase Kynurenine 3-hydroxylase KMO
Observed MW:	56 kDa
Uniprot ID :	O15229

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:100 - 1:400; IF: 1:10 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded mouse heart tissue slide using FNab04668(KMO Antibody) at dilution of 1:200



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab04668(KMO antibody) at dilution of 1:1000