

LETM2 antibody

Product Information

Catalog No.:	FNab04756
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	leucine zipper-EF-hand containing transmembrane protein 2
Synonyms:	LETM1 domain-containing protein LETM2, mitochondrial LETM1 and EF-hand domain-containing protein 2 Leucine zipper-EF-hand-containing transmembrane protein 1-like LETM2
Observed MW:	None
Uniprot ID :	Q2VYF4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue using FNab04756(LETM2 Antibody) at dilution of 1:50