

## LSM7 antibody

### Product Information

Catalog No.:	FNab04877
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds specifically to the 3'-terminal U-tract of U6 snRNA and is probably a component of the spliceosome.

### Immunogen information


Immunogen:	LSM7 homolog, U6 small nuclear RNA associated
Synonyms:	U6 snRNA-associated Sm-like protein LSm7 LSM7
Observed MW:	12 kDa
Uniprot ID :	Q9UK45

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab04877(LSM7 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04877(LSM7 antibody) at dilution of 1:1000