

## LTBP1 antibody

### Product Information

Catalog No.:	FNab04885
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be involved in the assembly, secretion and targeting of TGFB1 to sites at which it is stored and/or activated. May play critical roles in controlling and directing the activity of TGFB1. May have a structural role in the extra cellular matrix(ECM).

### Immunogen information


Immunogen:	latent transforming growth factor beta binding protein 1
Synonyms:	Latent-transforming growth factor beta-binding protein 1 (LTBP-1) Transforming growth factor beta-1-binding protein 1 (TGF-beta1-BP-1) LTBP1
Observed MW:	150 kDa
Uniprot ID :	Q14766

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:50-1:500; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Fnab04885( LTBP1 antibody) at dilution of 1:200



human plasma tissue were subjected to SDS PAGE followed by western blot with FNab04885( LTBP1 Antibody) at dilution of 1:2000