

MFAP5 antibody

Product Information

Catalog No.:	FNab05147
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in hematopoiesis. In the cardiovascular system, could regulate growth factors or participate in cell signaling in maintaining large vessel integrity(By similarity). Component of the elastin-associated microfibrils(PubMed:8557636).

Immunogen information


Immunogen:	microfibrillar associated protein 5
Synonyms:	Microfibrillar-associated protein 5 (MFAP-5) MP25 Microfibril-associated glycoprotein 2 (MAGP-2) MFAP5 MAGP2
Observed MW:	20 kDa
Uniprot ID :	Q13361

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:3000; IHC: 1:200-1:1000
Image:	



Immunohistochemistry of paraffin-embedded human placenta using FNab05147(MFAP5 antibody) at dilution of 1:100



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab05147(MFAP5 antibody) at dilution of 1:600