

MIA antibody

Product Information

Catalog No.:	FNab05173
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.

Immunogen information

Immunogen:	melanoma inhibitory activity
Synonyms:	Melanoma-derived growth regulatory protein Melanoma inhibitory activity protein MIA
Observed MW:	None
Uniprot ID :	Q16674

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skin cancer slide using FNab05173(MIA Antibody) at dilution of 1:50