

MRPL47 antibody

Product Information

Catalog No.:	FNab05335
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	mitochondrial ribosomal protein L47
Synonyms:	Large ribosomal subunit protein uL29m 39S ribosomal protein L47, mitochondrial (L47mt, MRP-L47) Nasopharyngeal carcinoma metastasis-related protein 1 MRPL47 NCM1
Observed MW:	29 kDa
Uniprot ID :	Q9HD33

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon tissue slide using FNab05335(MRPL47 Antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab05335(MRPL47 Antibody) at dilution of 1:300